Processing nanopore reads with dada2

Preprocessing reads

Preparing reads

Dorado doesn't support demultiplexing of dual indexes on 5' and 3' ends. Additionally, ligated library reads can be either orientation. Most straightforward approach to solve this, is utilising flanking adapter sequences and cutadapt. Another alternative is to demultiplex index pairs in forward and reverse orientation, then build scripts to reverse complement and merge reverse reads.

Extracting forward reads to fastq file can be performed with following command

cutadapt -g t2f...rc(t2r) -O 13 -trimmed-only -m 1300 -M 1650 -o forward_out.fastq.gz raw_reads.fastq.gz

Extracting reverse reads

cutadapt -g t2r...rc(t2f) -O 13 -m 1300 -M 1650 –trimmed-only -o reverse_out.fastq.gz raw_reads.fastq.gz

Note! O, e, m and M parameters can be used to reduce chances of misaligned matches

Next step is to reverse complement reverse read file and join all reads to single file. This will make life easier later in the workflow.

seqkit seq -rp -seq-type DNA reverse_out.fastq.gz | gzip >reverse_comp.fastq.gz

Final step is to merge two files together

cat forward_out.fastq.gz reverse_comp.fastq.gz >nanopore_reads.fastq.gz

Demultiplexing reads

Prepare list of barcodes as a fasta file

Use cutadapt to demux fastq file. In example output files are written to demuxed subdirectory cutadapt -e 0 -O 12 -g file:barcodes.fasta -o "demuxed/{name}.fastq.gz" input.fastq.gz

Trimming PCR amplification primers

You can use cutadapt and bash scripts to trim forward and revers PCR primers from demultiplexed sequence files.

Importing set1

Loading libraries

```
library(dada2);packageVersion("dada2")
[1] '1.30.0'
library(knitr);packageVersion("knitr")
[1] '1.45'
library(Biostrings);packageVersion("Biostrings")
[1] '2.70.2'
library(DECIPHER);packageVersion("DECIPHER")
[1] '2.30.0'
library(phyloseq);packageVersion("phyloseq")
[1] '1.46.0'
library(tidyverse);packageVersion("tidyverse")
[1] '2.0.0'
library(kableExtra);packageVersion("kableExtra")
[1] '1.4.0'
library(mia);packageVersion("mia")
[1] '1.10.0'
# Path variables
" rath Variables
path <- "data/processed/set2"
training <- "~/feature_classifiers/SILVA_SSU_r138_2019.RData"
meta_file <- "data/set2_meta.tsv"
exportloc <- "results_set2/"
# Variables: truncation length, phix (Illumina)
truncation <- 1400
#Creater results_directory</pre>
#Creates results directory
dir.create(exportloc)
#metadata
metadata <- data.frame(read_tsv(meta_file))</pre>
```


```
#List files inside directory
list.files(path)
```


```
[1] "barcode07_trimmed_all.fastq.gz" "barcode08_trimmed_all.fastq.gz"
[3] "barcode09_trimmed_all.fastq.gz" "barcode10_trimmed_all.fastq.gz"
[5] "barcode11_trimmed_all.fastq.gz" "barcode12_trimmed_all.fastq.gz"
[7] "filtered"


# Forward fastq filenames have format: SAMPLENAME_R1_001.fastq
fnFs <- sort(list.files(path, pattern="_trimmed_all.fastq.gz", full.names = TRUE))
# Extract sample names, assuming filenames have format: SAMPLENAME_XXX.fastq
sample.names <- sapply(strsplit(basename(fnFs), "_"), `[`, 1)
```

Plotting sequence quality profile for sample pairs

```
# Base quality plot
prI <- plotQualityProfile(fnFs[1:2])
prI</pre>
```







Denoising

Filter sequence data

Filtering reads (maxEE ≈ 1 error/200 bp sequence is good starting point for nanopore)

Dereplicating sequences

Step condenses identical sequences saving computational time

```
derepFs <- derepFastq(filtFs, verbose = TRUE)</pre>
```


Learning error rates

```
# Forward read error rate
errF <- learnErrors(derepFs, multithread = TRUE)</pre>
```

266337400 total bases in 190241 reads from 1 samples will be used for learning the error rates.

Plotting error rate

Plotting error rate profile for forward reads
plotErrors(errF, nominalQ = TRUE)



Denoising

```
dadaFs <- dada(derepFs, err = errF, multithread = TRUE)

Sample 1 - 190241 reads in 165100 unique sequences.
Sample 2 - 117666 reads in 105089 unique sequences.
Sample 3 - 180368 reads in 154701 unique sequences.
Sample 4 - 189903 reads in 152282 unique sequences.
Sample 5 - 170359 reads in 149023 unique sequences.
Sample 6 - 167216 reads in 147312 unique sequences.
```

Building asv table

```
seqtab <- makeSequenceTable(dadaFs)
# Dimensions of ASV table
dim(seqtab)</pre>
```

[1] 6 1273

Chimera removal

[1] 6 791

Summary

Table 1: Summary table

	Input	Filtered	DenoisedF	Nonchimeric	N:o of variants
PR07	246798	190241	175206	166449	153
PR08	152496	117666	102717	87434	150
PR09	234825	180368	175081	169869	132
PR10	244959	189903	184551	176725	120
PR11	220979	170359	154301	140981	152
PR12	212645	167216	145282	127029	165

Taxonomy assignment

IdTaxa from DECIPHER package

Create phyloseq object

Sequence data is stored as taxa_names. We will store sequences as refseq and create numbered variant names

Remove non-bacterial taxa

```
pseq <- subset_taxa(pseq, Kingdom != is.na(Kingdom))</pre>
```

Write results to files

Abundance table is transponed and written as tsv file

```
#variant names in rows
ASV_names <- taxa_names(pseq)
#sample names will be columns
ASV_counts <- t(otu_table(pseq))
ASVdf <- (data.frame(ASV_names,ASV_counts))
#write
write_tsv(ASVdf, paste0(exportloc,"asvs.tsv"))</pre>
```

Likewise taxonomy table is saved as tsv

```
#variant names in rows
ASV_names <- taxa_names(pseq)
#taxonomy ranks in columns
taxonomy <- (data.frame(ASV_names, tax_table(pseq)))
#write
write_tsv(taxonomy,paste0(exportloc,"taxonomy.tsv"))</pre>
```

Variant sequences are saved into fasta file