Processing nanopore reads

Marko Suokas

Preprocess reads

Dorado doesn't support demultiplexing of dual indexes on both 5' and 3' ends. Additionally, in ligated libraries reads can be either orientation. Our approach to demultiplex reads is using cutadapt. Index pairs are searched using linked adapters approach in forward and reverse orientation, then scripts are used to reverse complement reverse reads. Finally reads are merged

Please note that autocorrect often change double dashes in command examples.

Extracting forward reads to fastq file can be performed with following command

cutadapt -e 0 -O 12 -g file:~/scripts/barcodes.fasta –trimmed-only -m 1200 -o "fde-muxed/{name}.fastq.gz reads.fastq.gz

Command will extract barcodes defined in barcodes.fasta file and output matches into individual files in fdemuxed subdirectory. Minimum length is set in example to 1200 bp.

Extracting reverse reads using reverse complemented barcodes.fasta file

cutadapt -e 0 -O 12 -g file:~/scripts/rev_barcodes.fasta –trimmed-only -m 1200 -o "rde-muxed/{name}.fastq.gz reads.fastq.gz

Reads are demultiplexed into separate directory

Tip! O, e, m and M parameters can be used to reduce chances of misaligned matches

Next we use bash script that will process each reverse read file and reverse complement them using basic command

seqkit seq -rp -seq-type DNA -o reverse_comp.fastq.gz reverse_out.fastq.gz


Final step is merge. You can use simple bash script that will merge files with same base name from two separate directories using basic command

zcat forward_out.fastq.gz reverse_comp.fastq.gz >merged_reads.fastq.gz

Finally, you can use cutadapt and bash scripts to trim forward and reverse PCR primers from sequence reads.

Import set1 to R

Load libraries



Set variables

```
# Path variables
path <- "data/processed/set1"
training <- "~/feature_classifiers/SILVA_SSU_r138_2019.RData"
silva <- "~/feature_classifiers/silva_nr99_v138.1_train_set.fa.gz"
species <- "~/feature_classifiers/silva_species_assignment_v138.1.fa.gz"
meta_file <- "data/set1_meta.tsv"
exportloc <- "set1/"
# Variable truncation length
truncation <- 1400
#Creates results directory
dir.create(exportloc)
#metadata file to df
metadata <- read_tsv(meta_file, show_col_types = F)
metadata <- column_to_rownames(metadata, var = "Sampleid")</pre>
```

For project, we took advantage of computing power of CSC and imported already executed data objects. R code is unaltered. Execution is controlled by eval parameter in code chunk. RDS files also save resources and time when document is edited and checked.

```
#List files inside directory
list.files(path)


[1] "barcode01_trimmed.fastq.gz" "barcode02_trimmed.fastq.gz"
[3] "barcode03_trimmed.fastq.gz" "barcode04_trimmed.fastq.gz"
[5] "barcode05_trimmed.fastq.gz" "barcode06_trimmed.fastq.gz"
[7] "results"

# Forward fastq filenames have format: SAMPLENAME_R1_001.fastq
fnFs <- sort(list.files(path, pattern="_trimmed_all.fastq.gz", full.names = T))
# Extract sample names, assuming filenames have format: SAMPLENAME_XXX.fastq
sample.names <- sapply(strsplit(basename(fnFs), "_"), `[`, 1)</pre>
```

Plot sequence quality profile for samples

```
# Base quality plot
prsetI <- plotQualityProfile(fnFs[1:6])
prsetI

prsetI <- readRDS("rds/set1_rds/prsetI.rds")</pre>
```



Filter sequence data

Filtering reads (maxEE ≈ 1 error/200 bp sequence should be good starting point for this amplicon)


Learn error rates

```
# Forward read error rate
errF <- learnErrors(filtFs, multithread = T)

errF <- readRDS("rds/set1_rds/errF.rds")</pre>
```

Plot error rates

Plotting error rate profile for forward reads
plotErrors(errF, nominalQ = T)



Denoise

```
dadaFs <- dada(derepFs, err = errF, multithread = T)

dadaFs <- readRDS("rds/set1_rds/dadaFs.rds")</pre>
```

Build asv table

Dimensions tell us number of samples and variants

```
seqtab <- makeSequenceTable(dadaFs)
# Dimensions of ASV table
dim(seqtab)</pre>
```

[1] 6 72

Chimera removal

[1] 6 63

Summary

Table 1: Summary table

	Input	Filtered	DenoisedF	Nonchimeric	N:o of variants
barcode01	277408	225976	225599	223766	19
barcode02	277991	223708	223395	223395	12
barcode03	282347	221198	220613	220613	11
barcode04	261566	215868	215550	210277	8
barcode05	273330	203550	203212	203212	16
barcode06	261321	204075	203814	202829	8

Taxonomy assignment

Taxonomy classification against Silva 138.1 including species information.

```
taxonomy <- assignTaxonomy(seqtab.nochim, silva, multithread=3)
taxonomy <- addSpecies(taxonomy, species)
saveRDS(taxonomy, "rds/set1_rds/taxonomy.rds")

taxonomy <- readRDS("rds/set1_rds/taxonomy.rds")</pre>
```

Create TSE object

```
class: TreeSummarizedExperiment
dim: 63 6
metadata(0):
assays(1): counts
rownames(63): ASV1 ASV2 ... ASV62 ASV63
rowData names(7): Kingdom Phylum ... Genus Species
colnames(6): barcode01 barcode02 ... barcode05 barcode06
colData names(1): Name
reducedDimNames(0):
mainExpName: NULL
altExpNames(0):
rowLinks: NULL
rowTree: NULL
collinks: NULL
collree: NULL
collree: NULL
referenceSeq: a DNAStringSet (63 sequences)
```

Write results to files

Abundance table into tsv file

```
#sample names will be columns
ASVdf <- (data.frame(ASV_names,assays(tse_dada)$counts))
#write
write_tsv(ASVdf, paste0(exportloc,"asv_dada.tsv"))</pre>
```

Taxonomy table into tsv file

```
taxdf <- data.frame(ASV_names, rowData(tse_dada))
#write
write_tsv(taxdf, paste0(exportloc,"taxonomy_dada.tsv"))</pre>
```

Variant sequences into fasta file

Metadata into tsv file

```
metadf <- metadata %>% rownames_to_column(var = "Sampleid")
#write
write_tsv(metadf, paste0(exportloc,"metadata_dada.tsv"))
```

Add phylotree and save object

```
tree <- read.tree("set1/tree.nwk")
rowTree(tse_dada) <- tree
saveRDS(tse_dada, "set1/tse_dada.rds")</pre>
```

Vsearch@97%

Data has been processed in qiime, except taxonomic classification

```
class: TreeSummarizedExperiment
dim: 985 6
metadata(0):
assays(1): counts
rownames(985): ASV1 ASV2 ... ASV984 ASV985
rowData names(7): Kingdom Phylum ... Genus Species
colnames(6): barcode01 barcode02 ... barcode05 barcode06
colData names(1): Name
reducedDimNames(0):
mainExpName: NULL
altExpNames(0):
rowLinks: NULL
rowTree: NULL
collinks: NULL
colTree: NULL
referenceSeq: a DNAStringSet (985 sequences)
```

Write vsearch97 object

Vsearch@99%

```
class: TreeSummarizedExperiment
dim: 9195 6
metadata(0):
assays(1): counts
rownames(9195): ASV1 ASV2 ... ASV9194 ASV9195
rowData names(7): Kingdom Phylum ... Genus Species
colnames(6): barcode01 barcode02 ... barcode05 barcode06
colData names(1): Name
reducedDimNames(0):
mainExpName: NULL
altExpNames(0):
rowLinks: NULL
rowTree: NULL
colLinks: NULL
colTree: NULL
referenceSeq: a DNAStringSet (9195 sequences)
```

Write vsearch99 object

Observations

The low bacterial diversity in these samples may explain why denoising yields good results for long 16S rRNA sequences. The error rate plot appears flawless for this data. However, it is noteworthy that all samples contain over 150,000 unique reads.

In contrast, vsearch clustering generated a significantly higher number of variants, exceeding 900 and 9,000, respectively. Lowest number of variants (37) was observed with emu.